



LIFESTYLE DAY CHAIR BY MELROSE

The Melrose Hawk comes standard with a lightweight titanium frame that is nice and compact with no under-frame.

It is custom built to the individuals specifications for seat width and length. It has an adjustable camber bar position that allows you to alter your centre of gravity and seat height.

This makes it suitable for new wheelchair users that are unsure what seating position will work best for them.



A castor with adjustable positioning is necessary when adjusting the squeeze and rake

HAWK

STANDARD SPECIFICATIONS

- Titanium frame
- Adjustable back height and angle
- Fixed or folding back
- Adjustable footplate height
- Adjustable seat height
- Adjustable camber bar 0° to 5°
- Adjustable balance point
- Castors: rollerblade
4", 5", 6" solid
4", 5", 6" wide solid
- Quick release 24" wheels spoked on alloy rims
- Tyres high pressure tread or smooth, 24 x 1 solid KIK
- High or low mounted push to lock or sissor brakes
- Black fabric sling and back with leather trim
- Leg strap and fabric skirt guards
- Five-year frame warranty

OPTIONS AND ACCESSORIES

- Aluminum frame
- Castors: 4", 5", 6" solid alloy
4", 5", 6" soft roll
6", 8" pneumatic
- Frog Leg forks
- Quad release axles
- Wheels: Spoked 24", 25", 26"
BMX wheels
Xcore wheels
Spinergy wheels
Spider wheels
- Pushrims: Rubber coated - Red or Black
Vinyl
Foam
Max Greppe
Natural fit
- Spoke guards - clear or black
- Push handles: Removeable
- Brakes: Hill climber brakes
Superlock brakes
Sissor brakes
D-lock brakes
Brake extension levers
- Anti-tips with plastic wheels or rollerblades
- Armrest: Tubular swing away
Desk or full length pad
- adjustable height, with/without skirt guards
- Adjustable back upholstery
- ABS skirt guards or fabric covered ABS skirt guards
- Bag under chair
- Cushion
- Fabric or ABS footplate cover
- Impact guards, castor covers

