



panther

MELROSE LIFESTYLE WHEELCHAIR

FIXED FRAME DAY CHAIR

The Panther is one of our popular lightweight fixed frame wheelchairs. It is compact and ultra light, making it ideal for the active person.

A user friendly chair, it is also the lowest priced custom made wheelchair in the Melrose range. It has an adjustable camber bar position allowing you to alter your centre of gravity and it can also be built with the Ergo Seat System.



Piranha with fabric skirt guards

PANTHER

STANDARD SPECIFICATIONS

- Aluminum frame
- Adjustable back height and angle
- Fixed or folding back
- Adjustable footplate height
- Fixed camber and customer specified fixed seat height
- Castors: rollerblade
 - 4", 5", 6" solid
 - 4", 5", 6" wide solid
- Quick release 24" wheels spoked on alloy rims
- Tyres high pressure tread or smooth, 24 x 1 solid KIK
- High or low mounted push to lock or sissor brakes
- Black fabric sling and back with leather trim
- Leg strap and fabric skirt guards
- Five-year frame warranty

OPTIONS AND ACCESSORIES

- Titanium frame
- Ergo Seating
- Castors: 4", 5", 6" solid alloy
 - 4", 5", 6" soft roll
 - 6", 8" pneumatic
- Frog Leg forks
- Quad release axles
- Wheels: Spoked 24", 25", 26"
 - BMX wheels
 - Xcore wheels
 - Spinergy wheels
 - Spider wheels
- Pushrims: Rubber coated - Red or Black
 - Vinyl
 - Foam
 - Max Greppe
 - Natural fit
- Spoke guards - clear or black
- Push handles: Removeable
- Brakes: Hill climber brakes
 - Superlock brakes
 - Sissor brakes
 - D-lock brakes
 - Brake extension levers
- Anti-tips with plastic wheels or rollerblades
- Armrest: Tubular swing away
 - Desk or full length pad
 - adjustable height, with/without skirt guards
- Adjustable back upholstery
- ABS skirt guards or fabric covered ABS skirt guards
- Bag under chair
- Cushion
- Fabric or ABS footplate cover
- Impact guards, castor covers
- Changeable camber bar 0°, 3° or 6°

horizontal
seat

angled
seat

ergo
seat

MELROSE ERGO SEATING SYSTEM

A horizontal seat provides good pressure distribution with a neutral pelvic position while an angled seat gives better trunk stability and better pushing.

The ergonomic seat combines both horizontal and angled seating to give all three benefits. The ergonomic seat is angled, except for the back section, which is horizontal. The backrest is vertical and includes a rear-angled section at the top to give lumbar support for the spine.

