



### **FIXED FRAME DAY CHAIR**

The Melrose Piranha has a lightweight aluminum frame that is extra compact, even more compact than our Extreme.

It is completely custom built to the individuals specifications for size, rake and camber. It has an adjustable camber bar position allowing you to alter your centre of gravity. It can also be built with the Ergo Seat System



*Fixed camber bar position Piranha with impact guards, aluminium skirt guards, and push to lock super retractable brakes.*

## PIRANHA

### STANDARD SPECIFICATIONS

- Aluminum frame
- Adjustable back height
- Fixed or folding back
- Adjustable footplate height
- Adjustable or fixed camber bar position
- Customer specified fixed seat height
- Castors: rollerblade  
4", 5", 6" solid  
4", 5", 6" wide solid
- Quick release 24" wheels spoked on alloy rims
- Tyres high pressure tread or smooth, 24 x 1 solid KIK
- High or low mounted push to lock or sissor brakes
- Black fabric sling and back with leather trim
- Leg strap and fabric skirt guards
- Five-year frame warranty

### OPTIONS AND ACCESSORIES

- Titanium frame
- Ergo Seating
- Castors: 4", 5", 6" solid alloy  
4", 5", 6" soft roll  
6", 8" pneumatic
- Frog Leg forks
- Quad release axles
- Wheels: Spoked 24", 25", 26"  
BMX wheels  
Xcore wheels  
Spinergy wheels  
Spider wheels
- Pushrims: Rubber coated - Red or Black  
Vinyl  
Foam  
Max Greppe  
Natural fit
- Spoke guards - clear or black
- Push handles: Removeable
- Brakes: Hill climber brakes  
Superlock brakes  
Sissor brakes  
D-lock brakes  
Brake extension levers
- Anti-tips with plastic wheels or rollerblades
- Armrest: Tubular swing away  
Desk or full length pad  
- adjustable height, with/without skirt guards
- Adjustable back upholstery
- ABS skirt guards or fabric covered ABS skirt guards
- Bag under chair
- Cushion
- Fabric or ABS footplate cover
- Impact guards, castor covers

horizontal  
seat

angled  
seat

ergo  
seat

**MELROSE ERGO SEATING SYSTEM**

*A horizontal seat provides good pressure distribution with a neutral pelvic position while an angled seat gives better trunk stability and better pushing.*

*The ergonomic seat combines both horizontal and angled seating to give all three benefits. The ergonomic seat is angled, except for the back section, which is horizontal. The backrest is vertical and includes a rear-angled section at the top to give lumbar support for the spine.*

