

Bowls Script Form

Fixed or Folding Frame Sports Wheelchair

90875

To Order Fill in the frame dimensions and features required. For variations describe in detail on a separate page. If you have any further questions contact us at the address at the bottom of the page.


Date _____	Delivery details _____
Name _____	_____
Address _____	_____
_____	_____
Phone _____	_____
Fax / Cell _____	Required by _____
Email _____	(Delivery is normally 4-6 weeks from order confirmation)
Disability _____	Grading _____
Weight _____	Height _____
Shoe Size _____	T-Shirt Size _____

MEASUREMENTS (Refer to diagram below)

A Seat width	F Frame length (rear of backpost to centre of castor pivot)	K Castor Offset - right	Q <input type="checkbox"/> Fixed Frame 90875X
B Seat depth	G Seat to floor - front	M Camber - left (0°-5°)	R <input type="checkbox"/> Folding Frame 90875
C Back height	H Seat to floor - rear	N Camber - right (0°-5°)	Back Style (fixed frame wheelchair only)
D Back angle°	I Balance Point	O Wheel offset - left	<input type="checkbox"/> Folding Back <input type="checkbox"/> Fixed Back
E Seat to footplate	J Castor Offset - left	P Wheel offset - right	

ADDITIONAL NOTES

(you can type here)




Bowls Script Form

Fixed or Folding Frame Sports Wheelchair

Name _____


FOOT REST HANGER ANGLE

- ☐ 70°
☐ 80° (std)
☐ 90°




FOOT REST WIDTH

Width for Pro T and V style foot rests only




IMPORTANT
Foot rest Width
 When measuring
 across feet you must
 add 2" (50mm)

FIXED FRONT FOOT REST




IMPORTANT
 The styles in this box are only available
 on the Fixed Frame wheelchair




- ☐ Straight Hangers
☐ Pro T Hangers
☐ V Hangers
- ☐ Double Bow style foot plate (std)
☐ D style foot plate (std)
☐ Angle adjustable foot plate **AD11520**
☐ 2 piece, flip up plastic foot plates **AD11524**
☐ 1 piece, flip up angle adjustable alloy foot plate, opens to the left **AD11550**
☐ 1 piece, flip up angle adjustable alloy foot plate, opens to the right (pictured) **AD11550**
- ☐ ABS on foot plate **AD11500**
☐ 2 piece, flip up angle adjustable alloy foot plates **AD11521**

DROP IN FOOT REST



- ☐ Straight Hangers
☐ Pro T Hangers
- ☐ Double Bow foot plate **AD11580**
☐ D style foot plate **AD11580**
☐ Angle adjustable foot plate **AD11520**
☐ ABS on foot plate **AD11500**

SWING AWAY FOOT REST



- ☐ Straight Hangers
☐ Pro T Hangers
☐ V Hangers
- ☐ 2 piece, flip up plastic foot plates **AD11560**
☐ 1 piece, flip up angle adjustable alloy foot plate, opens to the left **AD11555**
☐ 1 piece, flip up angle adjustable alloy foot plate, opens to the right (pictured) **AD11555**
- ☐ 2 piece, flip up angle adjustable alloy foot plates **AD11570**

Bowls Script Form

Fixed or Folding Frame Sports Wheelchair

Name _____

CASTORS☐ 8" (200 x 50) pneumatic caster with tube wrap around (std)**FORKS**☐ Standard forks (std)**PUSH HANDLES**☐ Nil☐ In back canes (std) (folding frame wheelchair only)☐ Foldable in back canes **AD14160** (folding or fixed frame wheelchair)☐ Removable on grab bar **AD14100** (fixed frame wheelchair only)☐ Dual clamp adj height **AD14120** (fixed frame wheelchair only)**BRAKES**☐ Nil☐ Scissor (std)☐ Superlock **AD14560** (push to lock with brake extensions)☐ Brake extensions - left **AD14550**☐ Push to lock (std)☐ Hill climbers **AD14530**☐ Brake extensions - right **AD14550****ANTI-TIPS**☐ Nil☐ With plastic wheels **AD15000**☐ Adjustable anti tips **AD15010**☐ With roller blades **AD15030**☐ Kick tube **AD15040****ARMRESTS**☐ Nil☐ T adjustable height with desk pad armrest **AD15260**☐ ABS skirt guard attached to full length pad armrest **AD15281**☐ T adjustable height with full length pad armrest **AD15261****UPHOLSTERY**☐ Sling (std)☐ Standard back (std)☐ No cushion☐ Bag under chair **AD16600**☐ Tension adjustable back **AD16010**☐ Cushion 1" **AD16201**☐ Netting under chair **AD16650**☐ Cushion 2" **AD16202**☐ Bag on back of chair **AD16680**☐ Cushion 3" **AD16203**☐ Wheel bag **AD16780****SKIRT GUARDS**☐ Nil☐ ABS drop in **AD16110**☐ Alloy fixed **AD16130**☐ Carbon Fibre folding **AD16160**☐ Fabric skirt guards **AD16100**☐ ABS fabric covered **AD16120**☐ Alloy folding **AD16140****STRAPS / GLOVES / IMPACT GUARDS**☐ Leg strap (std)☐ Gloves☐ Impact guards front hangers **AD16400**☐ Heel loop straps (std)(size & style)☐ Impact guards back bars **AD16420****WHEELS**☐ 20" BMX (std)☐ Other: _____**TYRES**☐ 20" x 1.75" Smooth tyres with tube wrap around (std)

Other: _____

PUSHRIMS☐ Anodised (std)☐ Rubber Coated **AD13850**☐ Foam Coated **AD13840**☐ Vinyl Coated **AD13830****PUSHRIMS - POSITION**☐ Standard☐ In close☐ Onto rim**ACCESSORIES**☐ Ball holder (tray) **AD17020**☐ Ball Grabber **AD17021**

Bowls Script Form

Folding Frame Sports Wheelchair

Name _____

COLOURS

☐ 27 Candy Blue



☐ 1 Royal Blue



☐ 72 Glass Blue



☐ 16 Candy Purple



☐ 48 Burgundy



☐ 83 Candy Raspberry



☐ 88 Fluro Pink



☐ 28 Candy Red



☐ 11 Flame



☐ 57 Super Hot Chilli



☐ 80 Candy Orange



☐ 71 Pure Orange



☐ 68 Crystal Yellow



☐ 73 Crystal Jade



☐ 89 Fluro Green



☐ 70 Candy Emerald



☐ 2 Hawthorn Green



☐ 17 Crystal Sapphire



ATTENTION

THIS IS A GUIDE TO COLOURS ONLY, how colours appear on your monitor or when printed out may vary.

Contact us if you require assistance or a colour sample.

This colour selection is for our aluminum frames only.

We do not recommend powder coating titanium.

Colour availability is subject to change without notice, contact us if you are concerned about the availability of your colour.

☐ 55 Silver Sparkle



☐ 12 Silver Orbiter



☐ 46 Chrome



☐ 42 Metallic Viper



☐ 30 Satin Black



☐ 61 High Gloss Black




☐ 62 Matt Black



☐ 63 Texture Matt Black



☐ 18 Clear



☐ **Bead Blast** (titanium only)



☐ **Custom Colour** (additional cost)

Description