

Panther Script Form

Fixed Frame Day Wheelchair 90200

To Order Fill in the frame dimensions and features required. For variations describe in detail on a separate page. If you have any further questions contact us at the address at the bottom of the page.

Date _____		Delivery details _____	
Name _____		_____	
Address _____		_____	
_____		Therapist details _____	
Phone _____		_____	
Fax / Cell _____		Required by _____	
Email _____		TRIAL CHAIR REF: _____	
Disability _____	Weight _____	Height _____	Shoe Size _____ T-Shirt Size _____

MEASUREMENTS (Refer to diagram below)


Frame Materials

☐ Aluminium (std) ☐ Titanium **90200T**

A Seat Width _____	F Frame Length _____	K Castor Pivot Centres _____	<input type="checkbox"/> Ergo Seat Option AD10220
B Seat Length _____	G Seat to Floor - Front _____	L Footplate Width _____	N Ergo Seat Length _____
C Back Height _____	H Seat to Floor - Rear _____	M Camber (0° - 5°) _____	O Ergo Back Height _____
D Back Angle° _____	I Balance Point _____	Front Hanger Taper Style	P Ergo Back Angle° _____
E Seat to Footplate _____	J Castor Pivot Position _____	<input type="checkbox"/> V Front <input type="checkbox"/> Pro T Front	

ADDITIONAL NOTES

(you can type here)



Panther Script Form

Fixed Frame Day Wheelchair

Name _____

BACK REST

- ☐
- Folding back (std)
- ☐
- Fixed back

FOOTPLATE

- ☐
- Standard
- ☐
- Angle adjustable
- AD11520**
- ☐
- ABS on footplate
- AD11500**
- ☐
- Fabric padded cover
- AD11510**

CASTORS

- | | | | |
|---|---|---|--|
| <input type="checkbox"/> Rollerblade (std) | <input type="checkbox"/> 5" solid (std) | <input type="checkbox"/> 6" solid (std) | <input type="checkbox"/> 8" solid AD11760 |
| <input type="checkbox"/> 4" solid (std) | <input type="checkbox"/> 5" solid alloy AD11740 | <input type="checkbox"/> 6" solid alloy AD11742 | <input type="checkbox"/> 8" pneumatic AD11780 |
| <input type="checkbox"/> 4" x 1.25" alloy AD11710 | <input type="checkbox"/> 5" x 1.5" solid AD11720 | <input type="checkbox"/> 6" x 1.5" solid AD11730 | |
| <input type="checkbox"/> 4" wide soft roll AD11806 | <input type="checkbox"/> 5" x 1.5" alloy AD11740 | <input type="checkbox"/> 6" x 1.5" alloy AD11742 | |
| | <input type="checkbox"/> 5" wide soft roll AD11810 | <input type="checkbox"/> 6" wide soft roll AD11820 | |
| | | <input type="checkbox"/> 6" pneumatic AD11750 | |

FORKS

- ☐
- Standard forks (std)
- ☐
- Ultralight bearing forks
- AD11910**
- ☐
- Frog Legs suspension
- AD11950**
- ☐
- Frog Legs Unitine
- AD11920**

PUSH HANDLES

- ☐
- Nil
- ☐
- Removable
- AD14100**
- ☐
- Height adj dual clamp
- AD14120**
- ☐
- Foldable in back canes
- AD14160**

BRAKES

- | | | | |
|--|---|---|--|
| <input type="checkbox"/> Nil | <input type="checkbox"/> Superlock <small>Push to Lock with Brake Extensions</small> AD14560 | <input type="checkbox"/> Scissor (std) | <input type="checkbox"/> Scissor Lock Composite AD14570 |
| <input type="checkbox"/> Push to lock (std) | <input type="checkbox"/> Brake extension left AD14550 | <input type="checkbox"/> Ergo Scissor AD14522 | <input type="checkbox"/> D's-Locks AD14600 |
| <input type="checkbox"/> Pull to lock AD14510 | <input type="checkbox"/> Brake extension right AD14550 | <input type="checkbox"/> Ergo Scissor w/spring AD14524 | <input type="checkbox"/> Hill climbers AD14530 |

ANTI-TIPS

- | | |
|--|--|
| <input type="checkbox"/> Nil | <input type="checkbox"/> Single centre mount anti-tip with plastic wheels AD15070 |
| <input type="checkbox"/> Pair of single anti-tips with plastic wheels AD15000 | <input type="checkbox"/> Single centre mount anti-tip with rollerblades AD15080 |
| <input type="checkbox"/> Pair of single anti-tips with rollerblades AD15030 | <input type="checkbox"/> Dual tube one piece swivel anti-tip with rollerblade AD15150 |

ARMRESTS

- | | |
|---|--|
| <input type="checkbox"/> Nil | <input type="checkbox"/> T adjustable height with desk pad armrest AD15260 |
| <input type="checkbox"/> ABS skirt guard attached to full length pad armrest AD15281 | <input type="checkbox"/> T adjustable height with full length pad armrest AD15261 |

UPHOLSTERY

- | | | | |
|--------------------------------------|---|--|---|
| <input type="checkbox"/> Sling (std) | <input type="checkbox"/> Standard back (std) | <input type="checkbox"/> No cushion | <input type="checkbox"/> Bag under chair AD16600 |
| | <input type="checkbox"/> Tension adjustable back AD16010 | <input type="checkbox"/> Cushion 1" AD16201 | <input type="checkbox"/> Bag on back of chair AD16680 |
| | | <input type="checkbox"/> Cushion 2" AD16202 | <input type="checkbox"/> Bag on right front hanger AD16620 |
| | | <input type="checkbox"/> Cushion 3" AD16203 | <input type="checkbox"/> Bag on front left hanger AD16620 |
| | | | <input type="checkbox"/> Wheel bag AD16780 |

SKIRT GUARDS

- | | | | |
|---|--|---|--|
| <input type="checkbox"/> Nil | <input type="checkbox"/> ABS drop in AD16110 | <input type="checkbox"/> Alloy fixed AD16130 | <input type="checkbox"/> Carbon Fibre fixed AD16161 |
| <input type="checkbox"/> Fabric skirt guards AD16100 | <input type="checkbox"/> ABS fabric covered AD16120 | <input type="checkbox"/> Alloy folding AD16140 | <input type="checkbox"/> Carbon Fibre folding AD16160 |

STRAPS, IMPACT GUARDS

- | | | | |
|--|---|--|--|
| <input type="checkbox"/> Leg strap (std) | <input type="checkbox"/> Impact guards front hangers AD16400 | <input type="checkbox"/> Impact guards grab bar AD16420 | Impact guards made only if suitable after fitting anti tips, push handles etc. |
| | | <input type="checkbox"/> Impact guards bottom bar AD16422 | |

Panther Script Form

Fixed Frame Day Wheelchair

Name _____

WHEELS**24"**

- ☐ Spoked Wheels (std)
- ☐ Spinergy Day Wheels (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**
- ☐ Spinergy LX Wheels (12 spoke) **AD13310** (Choose your spoke colour) **Bk R Bl W Y Gr Pk O**
- ☐ Spinergy LX Wheels with FlexRim (12 spoke - black only) **AD13315**
- ☐ X Core Wheels (5 spoke) **AD13108**
- ☐ Spider Wheels **AD13100**
- ☐ Topolino Wheels **AD13350**
- ☐ Dino Wheels **AD13360** (Choose your spoke colour) **Bk R Bl Sv Go Gr Pk O LB Pu**
- ☐ Glance Wheels **AD13400** (Wheel Type) _____
- ☐ BMX Wheels (std)

TYRES**24"**

- ☐ 24" x 1" HP Smooth (std)
- ☐ 24" x 1" HP Tread (std)
- ☐ 24" x 1" Solid (std)
- ☐ 24" x 1-3/8" Tread (std)
- ☐ 24" x 1-3/8" Solid (std)
- ☐ 24" x 1-3/8" Kenda HP Tread **AD13710**
- ☐ 24" x 1-3/8" Knobbly V-Trak **AD13715**
- ☐ 24" x 1" Kenda Ironcap Blue **AD13700**
- ☐ 24" x 1" Kenda Konstrictor **AD13720**
- ☐ 24" x 1" Schwalbe Right Run **AD13740**
- ☐ 24" x 1" Schwalbe Marathon Evo **AD13750**
- ☐ BMX Tread (std)

25"

- ☐ Spoked Wheels (std)
- ☐ Spinergy Day (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**
- ☐ Spinergy LX (12 spoke) **AD13310** (Choose your spoke colour) **Bk R Bl W Y Gr Pk O**
- ☐ Spinergy LX Wheels with FlexRim (12 spoke - black only) **AD13315**
- ☐ X Core Wheels (3 spoke) **AD13200**
- ☐ Topolino Wheels **AD13350**
- ☐ Dino Wheels **AD13360** (Choose your spoke colour) **Bk R Bl Sv Go Gr Pk O LB Pu**
- ☐ Glance Wheels **AD13400** (Wheel Type) _____
- ☐ BMX Wheels (std)

25"

- ☐ 25" x 1" HP Smooth (std)
- ☐ 25" x 1" Solid (std)
- ☐ 25" x 1" Kenda Ironcap Blue **AD13700**
- ☐ 25" x 1" Kenda Konstrictor **AD13720**
- ☐ 25" x 1" Schwalbe Right Run **AD13740**
- ☐ 25" x 1" Schwalbe Marathon Evo **AD13750**
- ☐ BMX Tread (std)

26"

- ☐ Spoked Wheels (std)
- ☐ Spinergy Day (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**
- ☐ Spinergy LX (12 spoke) **AD13310**
- ☐ BMX Wheels (std)

26"

- ☐ 26" x 1" HP Smooth (std)
- ☐ 26" x 1-3/8" Tread (std)
- ☐ 26" x 1" Solid (std)
- ☐ 26" x 1" Schwalbe Right Run **AD13740**
- ☐ 26" x 1" Schwalbe Marathon Evo **AD13750**
- ☐ BMX Tread (std)

☐ Spare Wheel(s) Type: _____

Tyres: _____

Pushrims: _____

Size: _____

AXLES

- ☐ Quick release axles (std)
- ☐ Quad release **AD12500**

PUSHRIMS

- ☐ Anodised (std)
- ☐ Anodised (1") **AD13855**
- ☐ Foam **AD13840**
- ☐ Natural Fit Lite **AD13860**
- ☐ Vinyl **AD13830**
- ☐ Vinyl (1") **AD13859**
- ☐ Surge **AD13864**
- ☐ Max Greppe **AD13870**
- ☐ Rubber coated (black) **AD13850**
- ☐ Rubber coated (1") **AD13857**
- ☐ Surge Lite **AD13864**
- ☐ Spinergy FlexRim - see Wheels

PUSHRIMS - POSITION

- ☐ Standard
- ☐ In close
- ☐ Onto rim

SPOKE GUARDS

- ☐ Nil
- ☐ Clear polycarb **AD13900**
- ☐ Black ABS **AD13901**

Panther Script Form

Fixed Frame Day Wheelchair

Name _____

COLOURS

☐ 27 Candy Blue☐ 1 Royal Blue☐ 72 Glass Blue☐ 16 Candy Purple☐ 48 Burgundy☐ 83 Candy Raspberry☐ 88 Fluro Pink☐ 28 Candy Red☐ 11 Flame☐ 57 Super Hot Chilli☐ 80 Candy Orange☐ 71 Pure Orange☐ 68 Crystal Yellow☐ 73 Crystal Jade☐ 89 Fluro Green☐ 70 Candy Emerald☐ 2 Hawthorn Green☐ 17 Crystal Sapphire

ATTENTION


THIS IS A GUIDE TO COLOURS ONLY, how colours appear on your monitor or when printed out may vary.

Contact us if you require assistance or a colour sample.

This colour selection is for our aluminum frames only.

We do not recommend powder coating titanium.

Colour availability is subject to change without notice, contact us if you are concerned about the availability of your colour.

☐ 55 Silver Sparkle☐ 12 Silver Orbiter☐ 46 Chrome☐ 42 Metallic Viper☐ 30 Satin Black☐ 61 High Gloss Black☐ 62 Matt Black☐ 63 Texture Matt Black☐ 18 Clear☐ **Bead Blast** (titanium only)☐ **Custom Colour** (additional cost)

Description