



ORION
melrose basketball wheelchair

Do you need a light-weight super responsive frame designed to increase first push maneuverability?

Built from aluminium or titanium your computer aided designed frame will offer maximum rigidity. This directs all kinetic pushing forces directly to the wheels with no energy lost in frame flexibility.

We custom build each wheelchair to ensure a perfect fit, and to make sure that the athlete gets the very best out of their equipment.

We have some of the very best fabricators in the country producing our wheelchairs, and only the best components are used, if they don't meet our high standards then we make them ourselves.

The Orion Basketball Wheelchair, custom made just for you.



ORION

STANDARD SPECIFICATIONS

- Lightweight Aluminium Frame
- Complete custom fit - built to individual specifications for Size, Rake and Camber
- Camber 5° - 20°
- Quick Release 24", 25" or 26" Spoked Wheels on Alloy Rims
- Smooth High Pressure Tyres
- Rollerblade Castors
- Black Fabric Sling and Back
- 5 Year Frame Warranty

OPTIONS AND ACCESSORIES

- Super Lightweight Titanium Frame
- Tapered Seat Width and Back
- Double or Single Swivel Anti-Tips
- Skirt Guards: Fabric Skirt Guards
Alloy or Carbon Fibre Skirt Guards
- Adjustable Back Upholstery
- Straps and Bindings for Hips and Legs
- Cushion
- Smooth, High Pressure, Asymmetric Kenda Kaliente or Schwalbe Speedair Tyres
- Super light weight Spinerigy Wheels - 24", 25", 26" or 700C with colour option Carbon Cable Spokes
- Wheel Bag to accommodate four 24", 25" or 26" Wheels