



# Bowls

Built to comply with New Zealand Bowls regulations for using wheelchairs on the green.

Bowls holders can be mounted above as shown or below the chair. The low ride height puts you closer to the green and the wide castor positions allow your feet to clear the castor wheels. Wheel/frame offsets can also be built in to help place you over the bowl.

  
**Melrose**  
KIWI CONCEPT CHAIRS

What ever you can think of to improve your game, we can build for you.



## BOWLS

### STANDARD SPECIFICATIONS

---

- Aluminium frame
- Complete custom fit-built to individual specifications for size, rake and camber
- Fixed or folding frames as standard
- Adjustable back height in fixed frame
- Adjustable footplate height
- Quick release, 20" x 1.75" wheels
- 200mm x 50mm castors
- Skirt guards and foot strap
- Powdercoat finish

### OPTIONS AND ACCESSORIES

---

- Titanium frame
- Upholstery with adjustable back
- Rigid skirt guards
- Waist strap
- Back bar cover
- Tapered seat
- Rubber coated pushrims
- Spoke guards, ABS on footplate
- Anti tips
- Custom powder coat