

Airline Aisle Wheelchair

FOLDABLE ONBOARD AIRCRAFT

The Melrose Folding Aisle Chair is an onboard aisle chair that provides mobility about the cabin during long flights. It collapses down to be safely stowed into most aircraft lockers.

The frame is hand made from aircraft grade aluminium and is lightweight and rust resistant. Many of our components are made inhouse to ensure we maintain exact tolerances from the highest quality materials. Flame retardant soft trim materials can be used and can be certified to FAR23.853(a) and FAR25.853(a) App F Part I(a)1(ii) if required.

All four wheels are solid, these provide less rolling resistance for easier pushing and maneuverability up and down the aisles.

We can finish the chair in a large range of powder-coat colours and we will also have your airline or airport logos printed onto the seat back.

STANDARD SPECIFICATIONS

- Lightweight aluminum frame
- Dimensions:
 - Seat Height from floor 430mm
 - Seat Width 345mm
 - Seat Length 430mm
 - Seat Back Height from cushion 370mm
 - Maximum Chair Width 390mm
 - Weight 8.1kg
- Dimensions Folded:
 - Height 240mm
 - Width 390mm
 - Length 740mm
- Max user weight 120kg
- Solid roller-blade castors with 4" wide solid rear wheels
- Foot operated push to lock brake
- Restraint straps with composite clasp
- Five-year frame warranty

OPTIONS AND ACCESSORIES

- Fluro and Metallic colours
- Seatbelt style buckles on straps
- Airline or Airport logo on the seat back



Fits in most overhead lockers



Logo on Seat Back

